

QUICK
REFERENCE
GUIDE

AIRWAY MANAGEMENT PRODUCTS

Shiley™ Endotracheal Tubes



Medtronic
Further. Together

Shiley™ endotracheal tube catalog number charts

DESCRIPTION	QTY PER BOX	SIZE (MM)						
		1.0	1.5	2.0	2.5	3.0	3.5	4.0

VENTILATOR-ASSOCIATED PNEUMONIA (VAP) REDUCTION

TaperGuard™ evac oral endotracheal tube, Murphy eye	10							
TaperGuard™ evac oral endotracheal tube, Murphy eye, with pre-loaded stylet								
Evac oral endotracheal tube Seal Guard, Murphy eye	10							

BASIC CUFFED

TaperGuard™ with stylet, endotracheal tube, Murphy eye	10							
TaperGuard™ oral/nasal endotracheal tube, Murphy eye	10							
Oral/nasal endotracheal tube Seal Guard, Murphy eye	10							
Hi-Lo oral/nasal endotracheal tube, cuffed, intermediate, Murphy eye	10					86442	86443	86444
Hi-Lo oral/nasal endotracheal tube, cuffed, hooded Murphy tip with Murphy eye	10							
Lo-Pro oral/nasal endotracheal tube, cuffed, Murphy eye	10					86043		86045
Lo-Pro oral/nasal endotracheal tube cuffed, Magill tip	10					86003		86005

BASIC CUFFLESS

Oral/nasal endotracheal tube, cuffless, Murphy eye	10			86232	86233	86234	86235	86236
--	----	--	--	-------	-------	-------	-------	-------

SIZE (MM)											
4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0	9.5	10

			18860	18865	18870	18875	18880	18885	18890		
			18860S	18865S	18870S	18875S	18880S	18885S	18890S		
			110860	110865	110870	110875	110880	110885	110890		

	18750S	18755S	18760S	18765S	18770S	18775S	18780S	18785S	18790S		18710S
	18750	18755	18760	18765	18770	18775	18780	18785	18790		18710
	109850	109855	109860	109865	109870	109875	109880	109885	109890		109810
86445	86446	86447	86448	86449	86450	86451	86452	86453	86454	86455	86456
	86107	86108	86109	86110	86111	86112	86113	86114	86115		86117
86046	86047	86048	86049	86050	86051	86052	86053	86054	86055	86056	86057
	86007		86009	86010	86011	86012	86013	86014	86015	86016	86017

86237	86238	86239	86240	86241	86242						
-------	-------	-------	-------	-------	-------	--	--	--	--	--	--

Shiley™ endotracheal tube catalog number charts (cont'd.)

DESCRIPTION	QTY PER BOX	SIZE (MM)						
		1.0	1.5	2.0	2.5	3.0	3.5	4.0

RAE CUFFED AND CUFFLESS

Shiley™ oral RAE endotracheal tube with TaperGuard™ cuff	10							
Oral RAE endotracheal tube, cuffed, Murphy eye	10							86209
Oral RAE endotracheal tube, cuffless, Murphy eye	10							
Shiley™ nasal RAE endotracheal tube with TaperGuard™ cuff	10							
Nasal RAE endotracheal tube, cuffless, Murphy eye	10					86283	86284	86285

SPECIALTY CUFFED AND CUFFLESS

Oral/nasal endotracheal tube, cuffed reinforced, Murphy eye	5							
Endotracheal tube cuffless, reinforced, Magill tip	5			86541	86542	86543	86544	
Laser oral/nasal endotracheal tube, dual cuffed, Murphy eye	5							
Laser oral/nasal endotracheal tube, cuffless, Magill tip	5				86400	86401	86402	
Microlaryngeal oral/nasal endotracheal tube, cuffed, Murphy eye	10							86385

McGRATH™ MAC VIDEO LARYNGOSCOPE

DESCRIPTION	QTY PER BOX	SIZE (MM)	SKU NO.
Handle	1	N/A	301-000-000
McGRATH™ MAC blade size 1	50	1.0	350-072-000
McGRATH™ MAC blade size 2	50	2.0	350-017-000
McGRATH™ MAC blade size 3	50	3.0	350-005-000
McGRATH™ MAC blade size 4	50	4.0	350-013-000
McGRATH™ MAC X blade size 3	10	3.0	X3-003-000
250-minute battery pack	1	N/A	340-000-000

SIZE (MM)											
4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0	9.5	10

	76251	76255	76260	76265	76270	76275	76280	76285	76290		
86199											
86266	86267	86268	86269	86270	86271						
			96360	96365	96370	96375	96380				
86286	86287	86288	86289	86290	86291						

			86548	86549	86550	86551	86552	86553	86554	86555	
86545	86546	86547									
86397	86394	86395	86398								
	86387		86389								

Shiley™ endotracheal tube catalog number charts (cont'd.)

DESCRIPTION	QTY PER BOX	SIZE (MM)		
		24	26	28

ENDOBONCHIAL TUBES

Endobronchial tube, left	1			125028
Endobronchial tube, right	1			
Endobronchial tube, left with carinal hook	1			
Endobronchial tube, left with CPAP	1			125128
Right with CPAP system	1			
CPAP system for dual-lumen endotracheal tube (catalog no. 85006)	5			

ESOPHAGEAL/ENDOTRACHEAL DOUBLE-LUMEN AIRWAY

Use size 37 for patients 4' to 5'6" tall

Use size 41 for patients 5' tall and above

Esophageal/endotracheal double-lumen airway with X-ray line	4			
Esophageal/endotracheal double-lumen airway single line	4			
Esophageal/endotracheal double-lumen airway rollup kit line	4			
Esophageal/endotracheal double-lumen airway trainer, not for human use, mannequin demonstrations only	1			
Nasopharyngeal airway	10	8888-247015	8888-247023	8888-247031

SIZE (MM)						
30	32	34	35	37	39	41

	125032		125035	125037	125039	125041
			126035	126037	126039	126041
			125235	125237	125239	125241
	125132		125135	125137	125139	125141
			126135	126137	126139	126141

				5-18537		5-18541
				5-18237		5-18241
				5-18437		5-18441
						5-18141
8888-247049	8888-247056	8888-247064				



TaperGuard™ evac oral endotracheal tube, murphy eye

Incorporates Shiley™ evac technology which has been shown to reduce VAP.¹ Incorporates the TaperGuard™ cuff shown to reduce microaspiration by an average of 90% compared to the Shiley™ Hi-Lo cuff.²



TaperGuard™ endotracheal tube with stylet, murphy eye

The Shiley™ oral/nasal endotracheal tube with TaperGuard™ cuff and stylet gives you the tools you need in one sterile, easy-to-access pack — no preloading necessary.



Oral RAE endotracheal tube with TaperGuard™ cuff

Preformed curve design intended to reduce kinks and disconnects, maintain a patent airway, and position the circuits and connectors away from the surgical field. Features TaperGuard™ cuff technology.



Nasal RAE endotracheal tube with TaperGuard™ cuff

Preformed curve design positions circuits and connectors away from the surgical field without needing to bend or kink standard straight endotracheal tubes.




TaperGuard™ oral/nasal endotracheal tube, murphy eye

Unique taper-shaped cuff shown to reduce microaspiration by an average of 90% compared to the Shiley™ Hi-Lo cuff.²



Oral/nasal cuffless endotracheal tubes

Tubes feature distal tip reference marks and depth marks every centimeter to assist placement. Available in precut lengths. Smooth hooded tip with Murphy eye.



Microlaryngeal oral/nasal endotracheal tube

Provides greater access to surgical field during microlaryngeal surgical procedures. Standard tube internal and outside diameters, with longer tube length.



Endobronchial tube

For independent ventilation of either lung. Polyurethane endotracheal cuff is designed to reduce the risk of tearing during intubation. Available in left, right, and left with carinal hook styles.



Reinforced endotracheal tube

Maintains patent airway when patient's head and neck are extended or flexed, when the patient is prone, or when the patient will be moved.



Laser oral/nasal endotracheal tube

Designed for use with CO₂ and KTP laser surgery in close proximity to the airway. Flexible, stainless-steel body is airtight and laser resistant. Beams reflected from tube are defocused to help protect surrounding tissue.




Endobronchial tube with CPAP*

CPAP improves oxygenation during one-lung anesthesia. Adjustable valve provides clinically accurate pressure settings from 1–10 cm H₂O. Complete, ready-to-use system includes anesthesia bag.



Esophageal/endotracheal double-lumen airway

Unique design lets you establish a patent airway when placed in the trachea or esophagus.



Intubating stylet

Facilitates intubation of endotracheal tubes. 20 per box.

<u>Size</u>	<u>Catalog no.</u>
6 Fr	85863
10 Fr	85864
14 Fr	85865



Nasopharyngeal airway

Provides an open airway during suctioning and bronchoscopy. Soft material helps minimize patient coughing and gagging. Protects mucosa for repeated catheter passage.



McGRATH™ MAC enhanced direct laryngoscope

A single handheld device that combines the benefits of:

- Direct laryngoscopy, the familiar technique of a Macintosh laryngoscope
- Video laryngoscope video technology

*The CPAP valve is licensed by Dr. David Alfery.

1. Smulders K, van der Hoeven H, Weers-Pothoff I, Vandenbroucke-Grauls C. A randomized clinical trial of intermittent subglottic secretion drainage in patients receiving mechanical ventilation. *Chest*. 2002;121(3):858-862.
2. Shiley™ TaperGuard™ endotracheal tube and Shiley™ TaperGuard™ Evac endotracheal tube 510(k)#090352 benchtop testing.
3. Shiley™ endotracheal tube Seal Guard and Shiley™ endotracheal tube Seal Guard Evac 510(k)#082520 benchtop testing.

© 2016, 2017, 2019, 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 01/2020-16-AW-0068(4)-[WF# 2886722]

6135 Gunbarrel Avenue
Boulder, CO 80301

800.635.5267

[medtronic.com](https://www.medtronic.com)

Medtronic